

# Effects of anti-seizure compounds in a *DishBrain* system containing hyperactive glutamatergic neurons


Brad Watmuff<sup>1</sup>, Ben Rollo<sup>2</sup>, Koby Baranes<sup>3</sup>, Mark Kotter<sup>3</sup>, Brett J. Kagan<sup>1</sup>

1. Cortical Labs, Melbourne VIC AU. 2. Monash University Central Clinical School, Melbourne VIC AU. 3. University of Cambridge, UK.


*DishBrain*<sup>1</sup> is a novel *in vitro* model of neural information processing. Human neurons are presented structured information representing the game *Pong* and are able to interact with it through their activity.

Glutamatergic neuron hyperactivity is seen in pure glutamatergic cultures, and as a clinical feature is associated with a number of psychiatric and neurological conditions<sup>2</sup>.

We used *DishBrain* to determine whether modulation of gameplay performance via common anti-seizure drugs in a highly enriched glutamatergic culture is possible.




Measure neuronal responses to anti-seizure compounds



Generate dose-response relationships for anti-seizure compounds

Investigate effects of common anti-seizure medications on gameplay performance in the *DishBrain* platform


## Study design




## Anti-seizure compound effects on glutamatergic hyperactivity

Cultures set up in parallel to MEA chips were immunoreactive for neuronal and glutamatergic markers. MEA recordings revealed spikes corresponding with action potentials. Carbamazepine (200  $\mu\text{M}$ ) reduced mean and maximum firing rate significantly compared to time-matched vehicle control ( $p < 0.05$ ).


Induced NGN2 glutamatergic neurons




Neuronal firing (each dot = one action potential)



Anti-seizure compound effects on activity




Carbamazepine-activity dose-response curve




## *DishBrain* gameplay

Perampanel 1  $\mu\text{M}$  and carbamazepine 200  $\mu\text{M}$  improved hit/miss ratio and mean rally length in active gameplay sessions ( $p < 0.05$ ).

Anti-seizure compound effects on *Pong* hits



Anti-seizure compound effects on *Pong* rallies



## Key findings

Carbamazepine reduced the firing rate of highly enriched glutamatergic cultures. Dose-response relationships were generated with anti-seizure compounds in this assay. Perampanel and carbamazepine improved gameplay performance in the *DishBrain* system.

## Significance

We show varying performance in an *in vitro* neural processing model of glutamatergic hyperactivity using therapeutic interventions.

## Get in touch



## References

- Brett J. Kagan, Andy C. Kitchen, Nhi T. Tran, Bradyn J. Parker, Anjali Bhat, Ben Rollo, Adeel Razi, Karl J. Friston. In vitro neurons learn and exhibit sentience when embodied in a simulated game-world. bioRxiv 2021.12.02.471005; doi: <https://doi.org/10.1101/2021.12.02.471005>
- Miladinovic T, Nashed MG, Singh G. Overview of Glutamatergic Dysregulation in Central Pathologies. *Biomolecules*. 2015 Nov 11;5(4):3112-41. doi: 10.3390/biom5043112. PMID: 26569330; PMCID: PMC4693272.
- Pawlowski M, Ortmann D, Bertero A, Tavares JM, Pedersen RA, Vallier L, Kotter MRN. Inducible and Deterministic Forward Programming of Human Pluripotent Stem Cells into Neurons, Skeletal Myocytes, and Oligodendrocytes. *Stem Cell Reports*. 2017 Apr 11;8(4):803-812. doi: 10.1016/j.stemcr.2017.02.016. Epub 2017 Mar 23. PMID: 28344001; PMCID: PMC5390118.